CD9/CD9 Exosome ELISA Kit, Dog

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

Exosomes are membrane vesicles, which are secreted from many cells and have a diameter of about 30 nm ~ 200 nm. They are present in many bodily fluids, including saliva and blood^{1) ~ 5)}. It is also known that exosomes are secreted into the mammalian cell culture medium in vitro. Exosomes are enveloped with a lipid bilayer membrane, so like cells, membrane proteins exist at their surface, and they encase various molecular constituents, including proteins and microRNAs. When the target cells receive the exosomes, the proteins or microRNAs inside will fulfill their function, and in this way exosomes are thought to play a role of cell-to-cell communication⁶⁾. One of the exosomes' structural characteristics is the tetraspanin family which is located on their surface. Even in dogs, CD9 and CD63 are member molecules of this family, and known to be exosome surface markers^{7,8)}.

The conventional way to quantify exosomes is indirect. For example, it measures encapsulated protein abundance or performs Nanoparticle Tracking Analysis to determine the size distribution profile of small particles⁹⁾. The drawback for these methods is that they require exosome purification utilizing ultracentrifugation, for instance. The direct methods to quantitate exosomes in body fluids or culture supernatant is extremely limited, and there were no common methods available until now. The CD9/CD9 Exosome ELISA Kit is a Sandwich ELISA kit, which utilizes high-performance anti-dog CD9 antibodies. This product detects exosome markers, dog CD9 molecules that are located on the exosome surface in the body fluids or cell culture supernatant.

【 I – 2 】 Features

- Directly quantitate exosomes in human blood samples or cell culture supernatant.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Stabilized AZ-ACU (Canine urothelial carcinoma cells) derived exosomes are used as standard reagents.
- Normalization with Exosome Standard enable to relative quantitate each samples.
- Capture exosomes with solid phase anti-dog CD9 antibody, then detect using HRP conjugated anti-dog CD9 antibody.

【 I – 3 】 Kit Principle

This ELISA kit uses two-step sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-dog CD9 antibody.

First, samples were added onto the plate to capture exosomes by the anti-dog CD9 antibody. Next, HRP conjugated anti-dog CD9 antibody will be added to react exosome surface antigen, dog CD9. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-dog CD9 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	Exosome Standard	Dried product	4tubes* ¹
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X) * ²	25mL	1vial
5	HRP Conjugated Anti-dog CD9 Antibody (500X) * ³	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		2sheets

*¹ Sufficient to create 4 standard curves with n=2.

*² Crystals may precipitate in the Washing Buffer (10x) during refrigerated storage. Warm the buffer to dissolve it at 45°C before use.

*³ If the kit is not going to be used immediately, remove the labeled antibody from the kit and store it at -20°C.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μ L)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【Ⅱ】 Preparation of Reagents and Samples

【Ⅱ – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10 $\times$) to 10 folds with purified water.
e.g. For 1 plate, add 225mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【Ⅱ – 2】 Preparation of Exosome Standard solution

	Concentration (ng/mL)	Exosome Standard	Assay Buffer	Dilution factor
A	200	1vial	500 μ L	
B	100	250 μ L of A	250 μ L	2
C	50	250 μ L of B	250 μ L	2
D	25	250 μ L of C	250 μ L	2
E	12.5	250 μ L of D	250 μ L	2
F	6.25	250 μ L of E	250 μ L	2
G	3.13	250 μ L of F	250 μ L	2

- To prepare Solution A, add 500 μ L of Assay Buffer into 1 vial of Exosome Standard, and then mix well. To prepare Solution B, add 250 μ L of Assay Buffer into 250 μ L of Solution A, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 200 μ L for measurement, using 2 wells for each solution (n=2).

- Diluted Exosome Standard Solution(3.13~200 ng/mL) should be freshly prepared at each time before use.

【Ⅱ-3】 Preparation of antibody solution

- Dilute HRP conjugated anti-dog CD9 antibody (500x) to 500 folds using Assay Buffer.
e.g. For 1 plate, add 20μL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.
- * Diluted antibody solution should be freshly prepared at each time before use.

【Ⅱ-4】 Preparation of Samples (For Serum samples)

Serum is measured as a sample diluted 40-fold with Assay Buffer.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【Ⅱ-5】 Preparation of Samples (For cell culture medium supernatant)

The culture medium containing 10% fetal bovine serum (FBS) can be directly used to quantitate exosomes. Grow cells until it reach to the confluent state. Collect the medium supernatant, centrifuge at 2000 xg for 10min, and then used as a sample. If necessary, dilute with Assay Buffer.

【Ⅱ-6】 Sample Storage

Samples should be stored at 2~8℃ after preparation until measurement.

【Ⅲ】 Sample measurement procedure


1. Bring anti-dog CD9 antibody solid phased plate and the reagents to the room temperature.
2. Prepare dog CD9 Standard Beads solution by serial dilution. (from step 【Ⅱ-2】)
3. Add 100μL each of serial diluted Exosome Standard solution (3.13~200 ng/mL) or Sample solution into the well.

4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 2hrs for static reaction.
6. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100 μ L each of diluted HRP conjugated anti-dog CD9 antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake it in the plate shaker at 800 rpm for 30 sec.
9. Incubate at room temperature for 2hrs for static reaction.
10. Discard the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
11. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the average blank control subtracted absorbance value for each Exosome Standard concentration (y axis) against the Exosome Standard concentration (x axis).
15. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Determine the dog CD9 exosome concentration (Unit/mL) in the sample solution as the Exosome Standard (ng/mL) equivalent. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of dog CD9 exosome in the sample.

【IV】 Measurement example

【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against Exosome Standard concentration is drawn as shown in Figure 1. However, draw a new standard curve for each assay to calculate the concentration in the sample.

(Plotted is the value obtained by subtracting the blank absorbance from the absorbance of each standard protein concentration)

Exosome Standard (ng/mL)	Absorbance (450nm)		Average	Calculated concentration (U/mL)
	1	2		
0	0.115	0.103	0.109	
3.13	0.179	0.170	0.175	
6.25	0.213	0.219	0.216	
12.5	0.331	0.339	0.335	
25	0.546	0.581	0.564	
50	1.013	0.992	1.003	
100	1.767	1.834	1.801	
200	2.870	2.976	2.923	
Serum	0.302	0.315	0.309	11.5

Fig.1 Standard curve and measured values

【V】 Kit expiry date and storage

Expiry date : 6 months after the manufacturing date.

(The manufacturing date is indicated on the kit box label)

Storage : Refrigeration (2-8°C)

【Reference】

- 1) J. Skog, T. Wurdinger, S. van Rijn, D. H. Meijer, L. Gainche, M. Sena-Esteves, W. T. Jr. Curry, B. S. Carter, A. M. Krichevsky and X. O. Breakefield: *Nat Cell Biol.*, 10, 1470 (2008).
- 2) T. Pisitkun, R. F. Shen and M. A. Knepper: *Proc Natl Acad Sci USA.*, 101, 13368 (2004).
- 3) S. Runz, S. Keller, C. Rupp, A. Stoeck, Y. Issa, D. Koensgen, A. Mustea, J. Sehouli, G. Kristiansen and P. Altevogt: *Gynecol Oncol.*, 107, 563 (2007).
- 4) S. Keller, A. K. Konig, F. Marme, S. Runz, S. Wolterink, D. Koensgen, A. Mustea, J. Sehouli and P. Altevogt: *Cancer Lett.*, 278, 73 (2009).
- 5) C. Lasser, V. Seyed Ali, S. Gabrielsson, J. Lotvall and H. Valadi: *J Transl Med.*, 9, 9 (2011).
- 6) Y. Naito, Y. Yoshioka, Y. Yamamoto and T. Ochiya: *Cell Mol Life Sci.*, 74, 697 (2017).
- 7) E. J. Fish, K. J. Irizarry, P. DeInnocentes, et al.: *BMC Cancer.*, 18, 832 (2018).
- 8) M. A. Rojas, B. B. Rentzsch, J. Plendl, et. al.: *BMC Vet Res.*, 14, 179 (2018).
- 9) V. Filipe, A. Hawe and W. Jiskoot: *Pharm Res.*, **27**, 796 (2010).