# CDCP1 ELISA Kit, Human

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---

#### **【 I – 1 】 Background and Measurement Principal**

Prostate cancer (PCa) is one of the most frequent cancer types in men with 679,000 diagnoses and 220,000 deaths worldwide each year.

Most notably, several kallikrein-related peptidases, including KLK3 (prostate-specific antigen, PSA) have been extensively studied and found to promote prostate cancer progression. Recent findings also suggest a critical role for CDCP1 (CUB domain-containing protein 1) as the protease-activated transmembrane proteins in the development of advanced and aggressive castration-resistant prostate cancer (CRPC).<sup>1)</sup>

CDCP1 is overexpressed in a subset of CRPC. Notably, CDCP1 cooperates with the loss of the tumor suppressor gene PTEN to promote the emergence of metastatic prostate cancer. Moreover, we demonstrate that anti-CDCP1 immunoliposomes (anti-CDCP1 ILs) loaded with chemotherapy suppress prostate cancer growth when administered in combination with enzalutamide.<sup>2)</sup>

The effect of EVs isolated from metastatic PCa cells on mature osteoclast formation in the presence of RANKL and used functional siRNA screening to determine that CDCP1, which is a membrane protein located on EVs derived from PCa cells, promotes osteoclastogenesis. In addition, it was found that high levels of CDCP1 were detected in plasma-derived EVs from PCa patients with bone metastasis.<sup>3)</sup>

This product is a two-step sandwich ELISA kit to detect CDCP1 protein using a couple of highly specific monoclonal antibodies against the extracellular domain of human CDCP1.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

## 【 I – 2 】 Features

- Directly quantitate exosomal CDCP1 in human blood samples or cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Detect human CDCP1 by two-step sandwich method using immobilized anti-CDCP1 antibody and HRP conjugated anti-CDCP1 antibody.

## 【 I – 3 】 Kit Principle

This ELISA kit uses two-step sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-CDCP1 antibody.

First, the sample is added into the wells of the ELISA plate and allowed to react. Next after washing, HRP conjugated anti CDCP1 antibody is added to react.

After washing, substrate is added to the HRP conjugated anti CDCP1 antibody reacted with CDCP1. Finally, HRP color development is read with a plate reader to quantify human CDCP1 in the sample.

## 【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-CDCP1 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	CDCP1 Standard (20ng/mL)	150µL	1tube* <sup>1</sup>
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X) <sup>*2</sup>	25mL	1vial
5	HRP Conjugated Anti-CDCP1 Antibody (500X) <sup>*3</sup>	20µL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H <sub>2</sub> SO <sub>4</sub> )	6mL	1vial
8	Plate Seals		2sheet

\*<sup>1</sup> Sufficient to create 3 standard curves with n=2.

\*<sup>2</sup> Crystals may precipitate in the Washing Buffer (10x) during refrigerated storage. Warm the buffer to dissolve it at 45°C before use.

\*<sup>3</sup> If the kit is not going to be used immediately, remove the labeled antibody from the kit and store it at -20°C.

#### Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 µL)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

---

## 【II】 Preparation of Reagents and Samples

---

### 【II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to 10 folds with purified water.  
e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

### 【II – 2】 Preparation of Standard Protein solution

	Concentration (pg/mL)	CDCP1 Standard	Assay Buffer	Dilution factor
A	25,000			
B	2,500	50µL of A	450µL	10
C	1,250	250µL of B	250µL	2
D	625	250µL of C	250µL	2
E	312.5	250µL of D	250µL	2
F	156.3	250µL of E	250µL	2
G	78.2	250µL of F	250µL	2
H	39.1	250µL of G	250µL	2

- To prepare Solution B, add 450µL of Assay Buffer into 50µL of Standard Protein (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250µL of Assay Buffer

into 250µL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.

- Use 200µL for measurement, using 2 wells for each solution (n=2).
- Diluted CDCP1 Standard Solution(31.25~2,000pg/mL) should be freshly prepared at each time before use.

### **【II –3】 Preparation of antibody solution**

Dilute HRP conjugated anti-CDCP1 antibody (500x) to 500 folds using Assay Buffer.

e.g. For 1 plate, add 20µL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.

- \* Diluted antibody solution should be freshly prepared at each time before use.

### **【II –4】 Preparation of Samples**

Serum is measured as a sample diluted 4-fold with Assay Buffer.

Plasma is measured as a sample diluted 2-fold with Assay Buffer.

Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris.

Collect supernatants and assay.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

### **【II –5】 Sample Storage**

After sample preparation, store at 2-8°C until measurement.

---

## **【III】 Sample measurement procedure**

---

1. Bring anti-CDCP1 antibody solid phased plate and the reagents to the room temperature.
2. Prepare CDCP1 Standard solution by serial dilution. (from step 【II –2】 )
3. Add 100µL each of serial diluted CDCP1 Standard solution (31.25 ~ 2,000pg/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals.
5. Incubate at room temperature for 1 hour on a plate shaker set to 800 rpm.

6. Discard all the reaction solution, and then rinse each well with 300 $\mu$ L of Washing Buffer (from step 【II-1】 ). Repeat this step for 3 times.
7. Add 100 $\mu$ L each of diluted HRP conjugated anti CDCP1 antibody (from step 【II-3】 ) to the well.
8. Seal the microplate with Plate Seals.
9. Incubate at room temperature for 1 hour on a plate shaker set to 800 rpm.
10. Discard the reaction solution, and then rinse each well with 300 $\mu$ L of Washing Buffer (from step 【II-1】 ). Repeat this step for 3 times.
11. Add 100 $\mu$ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 $\mu$ L each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the absorbance value (y axis) against the Standard Protein concentration (x axis).
15. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of CDCP1 in the sample.


---

## 【IV】 Measurement example

---

### 【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against standard concentration is drawn as shown in Figure 1. However, draw a new standard curve for each assay to calculate the concentration in the sample.

Standard CDCP1 Solution (pg/mL)	Absorbance (450nm)		mean
	1	2	
0	0.056	0.056	0.056
31.25	0.106	0.105	0.106
62.5	0.155	0.149	0.152
125	0.252	0.262	0.257
250	0.454	0.445	0.450
500	0.858	0.815	0.837
1,000	1.660	1.508	1.584
2,000	2.724	2.535	2.630

Fig.1 Standard curve

---

### **【V】 Kit expiry date and storage**

---

Expiry date : indicated on the kit box

Storage : Refrigeration (2-8°C)

---

### **【Reference】**

---

- 1) H. Koistinen et al.: *IUBMB Life.*, **75**, 493 (2023).
- 2) A. Alajati et al.: *J Clin Invest.*, **130**, 2435 (2020).
- 3) F. Urabe et al.: *J Extracell Vesicles.*, **12**, 12312 (2023).