EphA2 ELISA Kit, Human

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EphA2 ELISA Kit, Human

【 I 】 About this kit

Erythropoietin-producing hepatocellular A2 (EphA2) is a member of the Eph tyrosine kinase receptor family that has been linked to various biological processes. In tumors, EphA2 overexpression is associated with noncanonical pathway activation, tumor progression, and a poor prognosis, which has emphasized its importance as a marker of malignancy.¹⁾

【 I – 1】 Features

- Consists of two monoclonal antibodies against extracellular domain of human EphA2 protein.
- Detects with high sensitivity Human EphA2 by sandwich method using a solid-phased anti-EphA2 antibody and an HRP-conjugated anti EphA2 antibody.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.

【 I – 2】 Kit Principle

This ELISA kit uses a two-step Sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti EphA2 antibody.

First, the sample is added into the wells of the ELISA plate and allowed to react.

Next after washing, HRP conjugated anti EphA2 antibody is added to react. After washing, substrate is added to the HRP conjugated anti EphA2 antibody reacted with EphA2.

Finally, HRP color development is read with a plate reader to quantify Human EphA2 in the sample.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

【 I – 3】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti EphA2 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	EphA2 Standard* ¹	200μL	1tube* ²
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X)* ³	25mL	1vial
5	HRP Conjugated Anti EphA2 Antibody (500X)* ⁴	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		1sheets

*¹ Purified recombinant human EphA2 (Gln25~Asn534) protein

*² Sufficient to create 4 standard curves with n=2.

*³ Crystals may precipitate in the Washing Buffer (10x) during refrigerated storage. Warm the buffer to dissolve it at 45°C before use.

*⁴ If the kit is not going to be used immediately, remove the labeled antibody from the kit and store it at -20°C.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μL)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【 II】 Preparation of Reagents and Samples

【 II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to 10 folds with purified water.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【II – 2】 Preparation of Standard Protein solution

	Concentration (pg/mL)	EphA2 Standard	Assay Buffer	Dilution factor
A	30000			
B	3000	50μL of A	450μL	10
C	1500	250μL of B	250μL	2
D	750	250μL of C	250μL	2
E	375	250μL of D	250μL	2
F	187.5	250μL of E	250μL	2
G	93.75	250μL of F	250μL	2
H	46.875	250μL of G	250μL	2

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of Human EphA2 Standard (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 100μL for measurement, using 2 wells for each solution (n=2).
- Diluted Human EphA2 Standard Solution(46.875~3000ng/mL) should be freshly prepared at each time before use.

【II – 3】 Preparation of antibody solution

- Dilute HRP conjugated anti-Human EphA2 antibody (500x) to 500 folds using Assay Buffer.

e.g. For 1 plate, add 20μL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.

* Diluted antibody solution should be freshly prepared at each time before use.

【II-4】 Preparation of Samples

In order to measure EphA2 protein dissolved in some buffers, the sample can be diluted with assay buffer appropriately. In case that EphA2 level in urine is analyzed, dilute the urine 10-fold with assay buffer and use as samples.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【II-5】 Sample Storage

After sample preparation, store at 2-8°C until measurement.

【III】 Sample measurement procedure

1. Bring anti-EphA2 antibody solid phased plate and the reagents to the room temperature.
2. Prepare Human EphA2 Standard solution by serial dilution. (from step 【II-2】)
3. Add 100μL each of serial diluted EphA2 Standard solution (46.875~3000pg/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals.
5. Incubate at room temperature for 1 hour on a plate shaker set to 800 rpm.
6. Discard all the reaction solution, and then rinse each well with 300μL of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100μL each of diluted HRP conjugated anti EPHA2 antibody (from step 【II-3】) to the well.
8. Seal the microplate with Plate Seals.
9. Incubate at room temperature for 1 hour on a plate shaker set to 800 rpm.
10. Discard the reaction solution, and then rinse each well with 300μL of Washing Buffer(from step 【II-1】). Repeat this step for 3 times.
11. Add 100μL of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50μL each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.


14. Create a standard curve by plotting the absorbance value (y axis) for each EphA2 Standard concentration against the EphA2 Standard concentration (x axis).
15. Determine the concentration of the target protein in the sample by interpolating absorbance values against the standard curve. Multiply the resulting value by the appropriate sample dilution factor, if used, to obtain the concentration of EPHA2 in the sample.

[IV] Measurement example

[IV- 1] Standard curve

As an example, the graph of absorbance (OD450) against Human Standard EphA2 concentration is drawn as shown in Figure 1.

However, draw a new standard curve for each assay to calculate the concentration in the sample.

EphA2 (pg/mL)	absorbance(450nm)		mean
	1	2	
0	0.050	0.061	0.056
46.875	0.095	0.099	0.097
93.75	0.146	0.149	0.148

187.5	0.243	0.251	0.247
375	0.465	0.445	0.455
750	0.853	0.842	0.848
1500	1.605	1.664	1.634
3000	2.839	2.737	2.788

Fig.1 Standard curve and measured values

【V】 Kit expiry date and storage

Expiry date : 6 months after the manufacturing date.

(The manufacturing date is indicated on the kit box label)

Storage : Refrigeration (2-8°C)

【Reference】

1) L. Toracchio, et al.: *Int J Mol Sci.*, **25**, 12191 (2024).