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
- (1) Alexandra C. Walls, Young-Jun Park, et al., Cell 180, 1-12 (2020)

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Created on December 21, 2023

Background

Angiotensin-converting enzyme 2 (ACE2), also called ACE homolog, is a membrane protein considerably homologous to ACE.

ACE2 was known as a factor in the renin-angiotensin system involved in blood pressure regulation.

Recently, it was clarified that the Spike glycoprotein of SARS-CoV-2 binds to ACE2 present on the cell membrane and then the virus is taken up into the cell when the infection occurs. Therefore, ACE2 is considered as a receptor for SARS-CoV-2 and this product (SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag) corresponds to the binding site of the Spike glycoprotein to the viral receptor ACE2⁽¹⁾.

As shown in the usage example, a highly sensitive binding assay can be performed by this product (HAK-SPD_bio-1) and human ACE2 His tag (HAK-ACE2_UL-1).

Product information

1. Protein Structure


Biotin-labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag (HAK-SPD_bio-1) contains Arg319-Phe541 of SARS-CoV-2 Spike glycoprotein tagged with Fc region of rabbit IgG1 and histidine at its C-terminus. The protein was expressed and secreted from HEK293 cells, followed by the purification using a Protein A column and then biotinylated.

2. Purity
> 95% (SDS-PAGE)

3. Composition
SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag 1mg/mL
0.1M PBS(pH7.2~7.4)
0.1% BSA

4. Storage:
Please avoid repeated freeze-thaw cycles.
The product must be stored at -70°C or lower upon receipt.
It is recommended to make small aliquots when using.

5. SDS-PAGE
500ng of the product before biotinylation was electrophoresed on 4~20% gradient gel and stained with CBB.

For research use only, not for diagnostic use.
Please read this manual carefully before use.

Usage example

1. Binding assay

Human ACE2 His tag was immobilized in a 96-well plate by adding at 0.1 μ g/100 μ L/well. Then 29-21000ng/mL of biotin labeled SARS-CoV-2 Spike glycoprotein Receptor-binding domain, Rabbit IgG1 Fc Tag was added to the wells. The binding was detected with HRP-labeled streptavidin. The linear range in the standard curve was 29-798ng/mL. Detail protocol for the binding assay is provided separately.