CD90/CD63 Exosome ELISA Kit, Human

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

Exosomes are membrane vesicles, which are secreted from many cells and have a diameter of about 30 nm ~ 200 nm. They are present in many bodily fluids, including saliva and blood ^{1) ~ 5)}. It is also known that exosomes are secreted into the mammalian cell culture medium in vitro. Exosomes are enveloped with a lipid bilayer membrane, so like cells, membrane proteins exist at their surface, and they encase various molecular constituents, including proteins and microRNAs. When the target cells receive the exosomes, the proteins or microRNAs inside will fulfill their function, and in this way exosomes are thought to play a role of cell-to-cell communication ⁶⁾. One of the exosomes' structural characteristics is the tetraspanin family which is located on their surface. CD9 and CD63 are member molecules of this family, and known to be exosome surface markers ⁷⁾.

Human mesenchymal stromal cells (hMSC) are multipotent cells with both regenerative and immunomodulatory activities making them an attractive tool for cellular therapy. Human mesenchymal stromal cells (hMSC) are multipotent cells with both regenerative and immunomodulatory activities making them an attractive tool for cellular therapy.^{8,9)}.

The conventional way to quantify exosomes is indirect. For example, it measures encapsulated protein abundance or performs Nanoparticle Tracking Analysis to determine the size distribution profile of small particles ¹⁰⁾. The drawback for these methods is that they require exosome purification utilizing ultracentrifugation, for instance. The direct methods to quantitate exosomes in body fluids or culture supernatant is extremely limited, and there were no common methods available until now.

This product is a two-step sandwich ELISA kit to quantitate the exosomes which are secreted from hMSC and positive in both CD63 which is a dominant tetraspanin in hMSC and CD90¹¹⁾

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

which is proposed as a marker for hMSC-derived exosomes utilizing high-performance anti-CD90 and anti-CD63 antibodies.

【 I – 2 】 Features

- Directly quantitate both CD90- and CD63-positve exosomes in cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.
- Utilize CD90/CD63 fusion protein (Standard Protein), instead of unstable/hard to store exosome itself, to implement stability and reproducibility.
- Normalization with a standard curve using CD90/CD63 fusion protein (Standard Protein) enable to relative quantitate each sample.
- Detect both CD90- and CD63-positve exosomes by two-step sandwich method using immobilized anti-CD90 antibody and HRP conjugated anti-CD63 antibody.

【 I – 3 】 Kit Principle

This ELISA kit uses two-step sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-CD90 antibody.

First, samples were added onto the plate to capture exosomes by the anti-CD90 antibody. Next, HRP conjugated anti-CD63 antibody will be added to react with CD63 on the surface of exosomes. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

| | Reagent | Volume | Quantity |
|---|---|-------------------------------|----------|
| 1 | Anti-CD90 Antibody Immobilized Plate | 96well (8well x 12 strips) | 1 plate |
| 2 | Standard Protein (CD90/CD63 Fusion Protein) | 200μL | 1tube* |
| 3 | Assay Buffer | 25mL | 1vial |
| 4 | Washing Buffer (10X) | 25mL | 1vial |
| 5 | HRP Conjugated Anti-CD63 Antibody (500X) | 20μL | 1tube |

| | | | |
|---|--|------|----------|
| 6 | Substrate Solution | 12mL | 1vial |
| 7 | Stop Solution (2N H ₂ SO ₄) | 6mL | 1vial |
| 8 | Plate Seals | | 2 sheets |

* Sufficient to create 4 standard curves with n=2.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 µL)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【II】 Preparation of Reagents and Samples

【II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to x10 with purified water.
e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【II – 2】 Preparation of CD9 Standard Beads solution

| | Concentration (pg/mL) | Standard Protein (CD90/CD63 Fusion Protein) | Assay Buffer | Dilution factor |
|---|--------------------------|---|--------------|--------------------|
| A | 8000 | | | |
| B | 800 | 50µL of A | 450µL | 10 |
| C | 400 | 250µL of B | 250µL | 2 |
| D | 200 | 250µL of C | 250µL | 2 |
| E | 100 | 250µL of D | 250µL | 2 |
| F | 50 | 250µL of E | 250µL | 2 |

| | | | | |
|---|------|------------|-------|---|
| G | 25 | 250μL of F | 250μL | 2 |
| H | 12.5 | 250μL of G | 250μL | 2 |

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of Standard Protein (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 100μL for measurement, using 2 wells for each solution (n=2).
- Diluted Standard Protein Solution (12.5~800pg/mL) should be freshly prepared at each time before use.

【II-3】 Preparation of antibody solution

Dilute HRP conjugated anti-CD63 antibody (500x) to 500 folds using Assay Buffer.

e.g. For 1 plate, add 20μL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.

* Diluted antibody solution should be freshly prepared at each time before use.

【II-4】 Preparation of Samples (For cell culture medium supernatant)

Centrifuge the cell culture supernatant at 2000 × g for 10 minutes to remove debris and use it as a sample. If dilution is necessary, dilute with assay buffer.


【II-5】 Sample Storage

Samples should be stored at 2~8°C after preparation until measurement.

【III】 Sample measurement procedure

1. Bring anti-CD90 antibody solid phased plate and the reagents to the room temperature.
2. Prepare Standard Protein solution by serial dilution. (from step 【II-2】)
3. Add 100μL each of serial diluted Standard Protein solution (12.5~800pg/mL) or Sample solution into the well.

4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.
6. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100 μ L each of diluted HRP conjugated anti-CD63 antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake in the plate shaker.
9. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.
10. Discard the antibody solution, and then rinse each well with 300 μ L of Washing Buffer. Repeat this step for 3 times.
11. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
13. Place into the plate reader, and read the absorbance of each microwell on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the average blank control subtracted absorbance value for each Standard Protein concentration (y axis) against the Standard Protein concentration (x axis) (Fig. 1).
15. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Determine the CD90-positive exosome concentration (Unit/mL) in the sample solution as the CD90/CD63 fusion protein (pg/mL) equivalent. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of both CD90- and CD63-positve exosomes in the sample.


| Standard Protein (CD90/CD63 Fusion Protein) (pg/mL) | Absorbance (450nm) | | Average |
|--|-----------------------|-------|---------|
| | 1 | 2 | |
| 0 | 0.060 | 0.056 | 0.058 |
| 12.5 | 0.097 | 0.097 | 0.097 |
| 25 | 0.149 | 0.146 | 0.148 |
| 50 | 0.248 | 0.241 | 0.245 |
| 100 | 0.447 | 0.446 | 0.447 |
| 200 | 0.859 | 0.851 | 0.855 |
| 400 | 1.649 | 1.611 | 1.630 |
| 800 | 2.984 | 2.977 | 2.981 |

Fig.1 Standard curve and measured values

【IV】 Sample measurement example

UE6E7T-3 cells, immortalized MSC, were seeded in a 12 well dish (1×10^4 cells/cm²) with 1mL POWEREDBY10 medium. The supernatant in each well was collected day by day, then centrifuged at 2,000 x g for five minutes to remove debris. The frozen stocks of the supernatants were thawed and directly measured with this kit.

Time course measurement of both CD90- and CD63-positive exosomes in immortalized MSC culture supernatant

【Reference】

- 1) J. Skog, T. Wurdinger, S. van Rijn, D. H. Meijer, L. Gainche, M. Sena-Esteves, W. T. Jr. Curry, B. S. Carter, A. M. Krichevsky and X. O. Breakefield: *Nat Cell Biol.*, **10**, 1470 (2008).
- 2) T. Pisitkun, R. F. Shen and M. A. Knepper: *Proc Natl Acad Sci USA.*, **101**, 13368 (2004).
- 3) S. Runz, S. Keller, C. Rupp, A. Stoeck, Y. Issa, D. Koensgen, A. Mustea, J. Sehouli, G. Kristiansen and P. Altevogt: *Gynecol Oncol.*, **107**, 563 (2007).
- 4) S. Keller, A. K. Konig, F. Marme, S. Runz, S. Wolterink, D. Koensgen, A. Mustea, J. Sehouli and P. Altevogt: *Cancer Lett.*, **278**, 73 (2009).
- 5) C. Lasser, V. Seyed AlikhS. Gabrielsson, J. Lotvall and H. Valadi: *J Transl Med.*, **9**, 9 (2011).

- 6) Y. Naito, Y. Yoshioka, Y. Yamamoto and T. Ochiya: *Cell Mol Life Sci.*, **74**, 697 (2017).
- 7) A. Zoraida and M. Yáñez-Mó: *Front Immunol.*, **5**, 442 (2014).
- 8) L. S. Meirelles, A. M. Fontes, D. T. Covas and A. I. Caplan: *Cytokine Growth Factor Rev.*, **20**, 419 (2009).
- 9) S. Rani, A. E. Ryan, M. D. Griffin and T. Ritter: *Mol Ther.*, **23**, 812 (2015).
- 10) V. Filipe, A. Hawe and W. Jiskoot: *Pharm Res.*, **27**, 796 (2010).
- 11) M. Théry et al.: *J Ext Vesicles*, **7**, 1535750 (2018).