EpCAM/CD9 Exosome ELISA Kit, Human

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

Human EpCAM protein is a transmembrane glycoprotein of 314 aa.

In the majority of cancer tissues, EpCAM is frequently overexpressed. Especially high abundant levels of EpCAM expression can be observed in carcinomas derived from colon, rectum, breast, lung and prostate.¹⁾ Besides affecting intercellular adhesion, EpCAM influences other important cellular functions relevant to tumor progression including cell proliferation and cancer stemness.²⁾

Exosomes are membrane vesicles which are secreted from many types of cells in the body and have a diameter of about 30nm~200nm. Recent reports demonstrate that most types of cancer tissue secrete large numbers of exosomes that carry molecular information about the parental tumor cells.³⁾ The information about exosomes derived from cancer cells could be an important tool for the clinical investigation.⁴⁾

However, as the analysis of exosome requires exosome purification processes using ultracentrifugation and so on, there were no rapid and convenient methods for direct measurement of exosome until now.

This product is a two-step sandwich ELISA, which utilizes high-performance anti-EpCAM and anti-CD9 monoclonal antibodies to directly and quantitatively detect EpCAM molecules on the surface of exosomes in human blood samples and cell culture supernatants.

【 I – 2 】 Features

- Directly quantitate exosomal EpCAM in human blood samples or cell culture supernatant.
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

- Utilize EpCAM/CD9 fusion protein (Standard Protein), instead of unstable/hard to store exosome itself, to implement stability and reproducibility.
- Normalization with a standard curve using EpCAM/CD9 fusion protein (Standard Protein) enable to relative quantitate each sample.
- Detect human EpCAM positive exosomes by two-step sandwich method using immobilized anti-CD9 antibody and HRP conjugated anti-EpCAM antibody.

【 I – 3 】 Kit Principle

This ELISA kit uses two-step sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-CD9 antibody.

First, samples were added onto the plate to capture exosomes by the anti-CD9 antibody. Next, HRP conjugated anti-EpCAM antibody will be added to react with EpCAM on the surface of exosomes. Finally, substrate will be added, then measure the coloring by the plate reader to quantitate sample exosomes.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-CD9 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	Standard Protein (EpCAM/CD9 Fusion Protein)	100μL	1tube*
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X)	25mL	1vial
5	HRP Conjugated Anti-EpCAM Antibody (500X)	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		2sheet

* Sufficient to create 2 standard curves with n=2.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μL)
- Multichannel micropipette
- Multichannel micropipette Reservoir

- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【II】 Preparation of Reagents and Samples

【II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to 10 folds with purified water.
e.g. For 1 plate, add 225mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【II – 2】 Preparation of Standard Protein solution

	Concentration (ng/mL)	Standard Protein	Assay Buffer	Dilution factor
A	200			
B	20	50μL of A	450μL	10
C	10	250μL of B	250μL	2
D	5	250μL of C	250μL	2
E	2.5	250μL of D	250μL	2
F	1.25	250μL of E	250μL	2
G	0.625	250μL of F	250μL	2
H	0.313	250μL of G	250μL	2

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of Standard Protein (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 200μL for measurement, using 2 wells for each solution (n=2).
- Diluted Standard Protein Solution(0.313~20 ng/mL) should be freshly prepared at each time before use.

【Ⅱ – 3】 Preparation of antibody solution

- Dilute HRP conjugated anti-EpCAM antibody (500x) to 500 folds using Assay Buffer.
e.g. For 1 plate, add 20μL of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.
- * Diluted antibody solution should be freshly prepared at each time before use.

【Ⅱ – 4】 Preparation of Samples (For Blood samples)

Serum or plasma is measured as a sample diluted 10-fold with Assay Buffer.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【Ⅱ – 5】 Preparation of Samples (For cell culture medium supernatant)

Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris.

Collect supernatants and assay.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【Ⅱ – 6】 Sample Storage

Store un-diluted sample at -20°C or below. Avoid repeated freeze-thaw cycles.

【Ⅲ】 Sample measurement procedure


1. Bring anti-CD9 antibody solid phased plate and the reagents to the room temperature.
2. Prepare Standard Protein solution by serial dilution. (from step 【Ⅱ – 2】)
3. Add 100μL each of serial diluted Standard Protein solution (0.156 ~ 10 ng/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.

6. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100 μ L each of diluted HRP conjugated anti-EpCAM antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake it in the plate shaker at 800 rpm for 30 sec.
9. Incubate at room temperature for 1 hour on a plat shaker set to 800 rpm.
10. Discard the reaction solution, and then rinse each well with 300 μ L of Washing Buffer(from step 【II-1】). Repeat this step for 3 times.
11. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the average blank control subtracted absorbance value for each Standard Protein concentration (y axis) against the Standard Protein concentration (x axis).
15. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Determine the EpCAM-positive exosome concentration (Unit/mL) in the sample solution as the EpCAM/CD9 fusion protein (ng/mL) equivalent. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of EpCAM positive exosome in the sample.

【IV】 Measurement example

【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against standard concentration is drawn as shown in Figure 1. However, draw a new standard curve for each assay to calculate the concentration in the sample.

(Plotted is the value obtained by subtracting the blank absorbance from the absorbance of each standard protein concentration)

Standard Protein (EpCAM/CD9 Fusion Protein) (ng/mL)	Absorbance (450nm)		Average	Calculated concentration (U/mL)
	1	2		
0	0.067	0.056	0.062	
0.313	0.099	0.098	0.099	
0.625	0.139	0.144	0.142	
1.25	0.222	0.224	0.223	
2.5	0.390	0.394	0.392	
5	0.745	0.746	0.746	
10	1.436	1.427	1.432	
20	2.716	2.689	2.703	
cell culture supernatants	0.979	1.085	1.032	7.07

Fig.1 Standard curve and measured values

【Reference】

- 1) PT. Went, A. Lugli, S. Meier, et al.: *Hum Pathol.*, **35**, 122 (2004).
- 2) K. Laura, W. Stefan, P. Klaus: *Cell Stress.*, **3**, 165(2019).
- 3) C. Théry, M. Ostrowski, E. Segura: *Nat Rev Immunol.*, **9**, 581 (2009).
- 4) AV. Vlassov, S. Magdaleno, R. Setterquist, et al.: *Biochim Biophys Acta.*, **1820**, 940 (2012).
- 5) C. Théry, S. Amigorena, G. Raposo, et al.: *Curr Protoc Cell Biol.*, **30**, Unit 3.22 (2006).