HER2 ELISA Kit, Dog

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【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

HER2 are frequently overexpressed in various human malignancies. Heterodimerization of these receptors results in the autophosphorylation of tyrosine residues within the cytoplasmic domain of the heterodimer and initiates signaling pathways leading to cellular proliferation and carcinogenesis.¹⁾

Lapatinib is an established therapeutic drug in human breast cancer harboring HER2 overexpression and/or gene amplification. It has been suggested that lapatinib (HER2 tyrosine kinase inhibitor) therapy combined with piroxicam in spontaneous canine muscle-invasive urothelial carcinoma is superior to piroxicam monotherapy in terms of clinical response, progression-free survival (PFS) and overall survival (OS).²⁾ As a result, it becomes important to detect the expression of HER2 in canine tumors as well.

In general, immunohistochemistry (IHC) tests for the protein expression in tissues are highly invasive and stressful for patients. In veterinary medicine, few reports describes quantitative measurement of HER2 protein in the canine body fluids although some IHC tests using anti-HER2 antibodies are reported.

This product is a two-step sandwich ELISA kit to detect HER2 protein in canine serum or urine using couple of highly specific monoclonal antibodies against canine HER2 which may be a potential biomarker in canine tumors.³⁾

【 I – 2 】 Features

- Detects high sensitivity dog HER2 by a two-step sandwich method using solid phase anti-dog HER2 antibody and HRP conjugated anti-dog HER2 antibody.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.

【 I – 3 】 Kit Principle

This ELISA kit uses two-step Sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-dog HER2 antibody.

First, the sample is added into the wells of the ELISA plate and allowed to react.

Next after washing, HRP conjugated anti-dog HER2 antibody is added to react. after washing, substrate is added to the HRP conjugated anti-dog HER2 antibody reacted with HER2.

Finally, HRP color development is read with a plate reader to quantify dog HER2 in the sample.

【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-dog HER2 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	Dog HER2 Standard	200μL	1tube*
3	Assay Buffer	25mL	1vial
4	Washing Buffer (10X)	25mL	1vial
5	HRP Conjugated Anti-dog HER2 Antibody (500X)	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H ₂ SO ₄)	6mL	1vial
8	Plate Seals		1sheets

* Sufficient to create 4 standard curves with n=2.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μL)
- Multichannel micropipette
- Multichannel micropipette Reservoir

- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

【II】 Preparation of Reagents and Samples

【II – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to 10 folds with purified water.
e.g. For 1 plate, add 225 mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【II – 2】 Preparation of Standard Protein solution

	Concentration (pg/mL)	Dog HER2 Standard	Assay Buffer	Dilution factor
A	30000			
B	3000	50μL of A	450μL	10
C	1500	250μL of B	250μL	2
D	750	250μL of C	250μL	2
E	375	250μL of D	250μL	2
F	187.5	250μL of E	250μL	2
G	93.8	250μL of F	250μL	2
H	46.9	250μL of G	250μL	2

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of Dog HER2 Standard (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 100μL for measurement, using 2 wells for each solution (n=2).
- Diluted Dog HER2 Standard Solution(46.9~3000 ng/mL) should be freshly prepared at each time before use.

【II –3】 Preparation of antibody solution

- Dilute HRP conjugated anti-dog HER2 antibody (500x) to 500 folds using Assay Buffer.
e.g. For 1 plate, add 20 μ L of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.
- * Diluted antibody solution should be freshly prepared at each time before use.

【II –4】 Preparation of Samples

Dilute serum 10-fold and urine 2-fold with assay buffer and use as samples. Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

【II –5】 Sample Storage

After sample preparation, store at 2-8°C until measurement.

【III】 Sample measurement procedure

1. Bring anti-dog HER2 antibody solid phased plate and the reagents to the room temperature.
2. Prepare Dog HER2 Standard solution by serial dilution. (from step 【II – 2】)
3. Add 100 μ L each of serial diluted Dog HER2 Standard solution (46.9~3000 pg/mL) or Sample solution into the well.
4. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
5. Incubate at room temperature for 2hrs for static reaction.
6. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
7. Add 100 μ L each of diluted HRP conjugated anti-dog HER2 antibody (from step 【II-3】) to the well.
8. Seal the plate, and then shake it in the plate shaker at 800 rpm for 30 sec.
9. Incubate at room temperature for 2hrs for static reaction.


10. Discard the reaction solution, and then rinse each well with 300 μ L of Washing Buffer(from step 【II-1】). Repeat this step for 3 times.
11. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
12. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
13. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
14. Create a standard curve by plotting the average blank control subtracted absorbance value for each Dog HER2 Standard concentration (y axis) against the Dog HER2 Standard concentration (x axis).
15. Determine the concentration of the target protein in the sample by interpolating the blank control subtracted absorbance values against the standard curve. Multiply the resulting value by the appropriate sample dilution factor, if used, to obtain the concentration of dog HER2 in the sample.

【IV】 Measurement example

【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against Dog Standard HER2 concentration is drawn as shown in Figure 1.

However, draw a new standard curve for each assay to calculate the concentration in the sample.

(Plotted is the value obtained by subtracting the blank absorbance from the absorbance of each standard protein concentration)

Dog HER2 (pg/mL)	Absorbance (450nm)		Average	Calculated value (pg/mL)
	1	2		
0	0.058	0.061	0.060	
46.9	0.101	0.099	0.100	
93.8	0.149	0.146	0.148	
187.5	0.246	0.254	0.250	
375	0.444	0.441	0.443	
750	0.885	0.816	0.851	
1500	1.593	1.563	1.578	
3000	2.967	2.900	2.934	
Serum sample	0.463	0.457	0.460	382.8

Fig.1 Standard curve and measured values

【Reference】

- 1) Y. Yarden, et al.: *Nat. Rev. Mol. Cell. Biol.* **2**, 127 (2001)
- 2) S. Maeda, et al.: *Sci Rep.*, **12**, 4 (2022)
- 3) Domestic patent pending.