# PD-L1 ELISA Kit, Human

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---

#### **【 I – 1 】 Background and Measurement Principal**

Tumor cells evade the immune surveillance by up-regulating surface expression of PD-L1, which interacts with PD-1 on T cells to elicit the immune checkpoint response<sup>1),2)</sup>. In recent years, anti-PD-1 antibody drugs and anti-PD-L1 antibody drugs that inhibit this signal transduction have been confirmed to have clinical efficacy and safety in various solid cancers by enhancing T cell activation<sup>3),4)</sup>.

However, therapeutic resistance after initial response has been increasingly observed, and mechanisms of resistance to anti-PD-L1 antibody therapies were mostly unknown.

Recently, it has been identified the unique secreted PD-L1 splicing variants lacking the transmembrane domain, from anti-PD-L1 therapies-resistant NSCLC patients<sup>5)</sup>.

It has been suggested that the secreted PD-L1 splicing variants induced the resistance to PD-L1 blockade therapy by acting as a “decoy”, and the levels of soluble PD-L1 in the exosome-free plasma or pleural effusion in patients detected with secreted PD-L1 splice variants were much higher than in healthy donors or patients without the secreted PD-L1 variants<sup>5)</sup>.

Increasing evidence suggests that the blood levels of soluble-PD-L1 might facilitate the prediction of clinicopathological characteristics, treatment response, and survival outcomes in the cancer patients<sup>6)</sup>.

This kit detects human PD-L1 by a one-step sandwich method using the high-performance antibody against human PD-L1, an immune checkpoint related molecule.

## 【 I – 2 】 Features

- Detects high sensitivity human PD-L1 with rapidly by a one-step sandwich method using solid phase anti-human-PD-L1 antibody and HRP conjugated anti-human-PD-L1 antibody. (Detection limit: 1pg / mL)
- No special equipment is required. Standard microplate reader capable of reading at 450nm will do the job.

## 【 I – 3 】 Kit Principle

This ELISA kit uses One-step Sandwich ELISA principle. The ELISA plate provided in this kit has been pre-coated with an anti-human PD-L1 antibody.

First, the sample is added into the wells of the ELISA plate with HRP conjugated anti-PD-L1 antibody and allowed to react. Next, after washing, substrate is added to the HRP conjugated anti-PD-L1 antibody reacted with PD-L1. Finally, HRP color development is read with a plate reader to quantify human PD-L1 in the sample.

## 【 I – 4 】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-PD-L1 Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
2	PD-L1 Standard	100μL	1tube*
3	Assay Buffer	15mL	1vial
4	Washing Buffer (10X)	15mL	1vial
5	HRP Conjugated Anti-PD-L1 Antibody (250X)	20μL	1tube
6	Substrate Solution	12mL	1vial
7	Stop Solution (2N H <sub>2</sub> SO <sub>4</sub> )	6mL	1vial
8	Plate Seals		1sheets

\* Sufficient to create 3 standard curves with n=2.

### Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000  $\mu\text{L}$ )
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)
- Plate washer

---

## 【 II 】 Preparation of Reagents and Samples

---

### 【 II – 1 】 Preparation of Washing Buffer(Prepare 2 wells for each assay)

- Dilute Washing Buffer (10 $\times$ ) to 10 folds with purified water.  
e.g. For 1 plate, add 135 mL of purified water to 15mL of Washing Buffer (10 x) and mix well.

### 【 II – 2 】 Preparation of Standard Protein solution

	Concentration (pg/mL)	PD-L1 Standard	Assay Buffer	Dilution factor
A	8000			
B	800	30 $\mu\text{L}$ of A	270 $\mu\text{L}$	10
C	400	150 $\mu\text{L}$ of B	150 $\mu\text{L}$	2
D	200	150 $\mu\text{L}$ of C	150 $\mu\text{L}$	2
E	100	150 $\mu\text{L}$ of D	150 $\mu\text{L}$	2
F	50	150 $\mu\text{L}$ of E	150 $\mu\text{L}$	2
G	25	150 $\mu\text{L}$ of F	150 $\mu\text{L}$	2
H	12.5	150 $\mu\text{L}$ of G	150 $\mu\text{L}$	2

• To prepare Solution B, add 270 $\mu\text{L}$  of Assay Buffer into 30 $\mu\text{L}$  of PD-L1 Standard (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 150 $\mu\text{L}$  of Assay Buffer into 150 $\mu\text{L}$  of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.

- Use 100 $\mu\text{L}$  for measurement, using 2 wells for each solution (n=2).

- Diluted PD-L1 Standard Solution(12.5~800 ng/mL) should be freshly prepared at each time before use.

### **【 II –3】 Preparation of antibody solution**

- Dilute HRP conjugated anti-PD-L1 antibody (250x) to 250 folds using Assay Buffer.  
e.g. For 1 plate, add 20 $\mu$ L of antibody (250x) into 5mL of Assay Buffer. Mix by inverting the tube.
- \* Diluted antibody solution should be freshly prepared at each time before use.

### **【 II –4】 Preparation of Samples (For Blood samples)**

Blood samples should be serum or plasma.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

### **【 II –5】 Preparation of Samples (For cell culture medium supernatant)**

Centrifuge cell culture media at 2,000 x g for 10 minutes to remove debris. Collect supernatants and assay.

Samples generating absorbance values greater than that of the highest standard should be further diluted using Assay Buffer and reanalyzed.

### **【 II –6】 Sample Storage**

Store un-diluted sample at -20 °C or below. Avoid repeated freeze-thaw cycles.

---

## **【 III】 Sample measurement procedure**

---

1. Bring anti-PD-L1 antibody solid phased plate and the reagents to the room temperature.
2. Prepare PD-L1 Standard solution by serial dilution. (from step 【 II – 2】 )
3. Add 50 $\mu$ L each of serial diluted PD-L1 Standard solution (12.5~800 pg/mL) or Sample solution into the well.

4. Furthermore, add 50 $\mu$ L each of diluted HRP conjugated anti-PD-L1 antibody (from step 【II-3】 ) to the well.
5. Seal the plate, and then shake in the plate shaker, and then shake it at 800 rpm for 30 sec.
6. Incubate at room temperature for 1 hour on a plate shaker set to 800 rpm, or for 2hours for static reaction.
7. Discard the reaction solution, and then rinse each well with 300 $\mu$ L of Washing Buffer(from step 【II-1】 ). Repeat this step for 3 times.
8. Add 100 $\mu$ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
9. Visually confirm the coloring, and then add 50 $\mu$ L each of Stop Solution.
10. Place into the Plate-reader, and read the absorbance of each well on a spectrophotometer at the wavelength of 450nm.
11. Create a standard curve by plotting the average blank control subtracted absorbance value for each PD-L1 Standard concentration (y axis) against the PD-L1 Standard concentration (x axis).
12. Determine the concentration of the target protein in the sample by interpolating the blank control subtracted absorbance values against the standard curve. Multiply the resulting value by the appropriate sample dilution factor, if used, to obtain the concentration of PD-L1 in the sample.

---


## 【IV】 Measurement example

---

### 【IV- 1】 Standard curve

As an example, the graph of absorbance (OD450) against standard PD-L1 concentration is drawn as shown in Figure 1.

However, draw a new standard curve for each assay to calculate the concentration in the sample.

PD-L1 (pg/mL)	Absorbance (450nm)		Average	Calculated value (pg/mL)
	1	2		
0	0.066	0.068	0.067	
12.5	0.106	0.113	0.110	
25	0.15	0.151	0.151	
50	0.239	0.240	0.240	
100	0.421	0.407	0.414	
200	0.789	0.783	0.786	
400	1.460	1.462	1.461	
800	2.690	2.595	2.643	
Serum sample	0.719	0.714	0.717	179.2

Fig.1 Standard curve and measured values

### 【IV-2】 Detection limit

The detection limit<sup>7)</sup> calculated from the standard deviation of the blank absorbance and the slope of the standard curve was 1 pg/mL.

### 【IV-3】 Recovery (Serum, Plasma)

Three concentrations of PD-L1 were spiked to evaluate signal recovery in the working range of the assay.

Table 1 Recovery of each sample

Sample	Recovery (%)	
	Average	Range
Serum	94%	92-99%
Plasma	91%	88-93%

---

### **【Reference】**

---

- 1) H. Dong, S. E. Strome, D. R. Salomao, H. Tamura, F. Hirano, D. B. Flies, P. C. Roche, J. Lu, G. Zhu, K. Tamada, V. A. Lennon, E. Celis, L. Chen: *Nat Med.*, **8**, 793 (2002).
- 2) L. Chen, X. Han: *J Clin Invest.*, **125**, 3384(2015).
- 3) S. L. Topalian, F. S. Hodi, J. R. Brahmer, et al.: *N Engl J Med.*, **366**, 2443 (2012).
- 4) J. R. Brahmer, S. S. Tykodi, L. Q. Chow, et al.: *N Engl J Med.*, **366**, 2455 (2012).
- 5) B. Gong, K. Kiyotani, S. Sakata, et al.: *J Exp Med.*, 216, 982 (2019).
- 6) X. Zhu, J. Lang: *Oncotarget*, 8, 97671 (2017).
- 7) Notification No.338 issued by the Director of Evaluation and Licensing Division, Pharmaceutical and Food Safety Bureau, Pharmaceutical and Medical Safety Bureau, Ministry of Health, Labour and Welfare, September 28, 1997.