SARS-CoV-2 Spike Glycoprotein ELISA Kit

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- (1) Alexandra C. Walls, Young-Jun Park, et al., Cell 180, 1-12 (2020)

SARS-CoV-2 Spike Glycoprotein ELISA Kit

【 I 】 About this kit

【 I – 1 】 Background and Measurement Principal

The pandemic caused by the new coronavirus infection (Covid-19), which spread from China to the world in 2019, have infected 197.19 million people and killed 4.21 million people by the end of July 2021 in the worldwide. SARS-CoV-2, which was identified as the causative virus of Covid-19, is a type of coronavirus that has an RNA with a total length of 30,000 base pairs as a gene. The spike protein encoded by the S gene in the viral genome is a membrane protein expressed on the surface of the virus. The spike protein is composed of two domains, S1 and S2.

It has been clarified that the S1 domain of the peplomer protein first binds to angiotensin converting enzyme 2 (ACE2) on the cell surface when SARS-CoV-2 infect cells, so that ACE2 is considered as a receptor for SARS-CoV-2. In addition, the structural analysis of the S1 domain has identified the narrow region which binds to ACE2 (receptor binding domain, RBD) (1).

This product can measure SARS-CoV-2 spike protein in a sample by the sandwich ELISA method using two types of neutralizing antibodies that are specific and have high affinity for the RBD.

【 I – 2 】 Features

- Measures SARS-CoV-2 spike protein in samples such as serum and plasma.
- No special equipment required, measurement is possible with a normal plate reader.
- SARS-CoV-2 spike protein can be measured sensitively and specifically by sandwich ELISA method using two types of antibodies against RBD in the SARS-CoV-2 S1 domain.

For research use only, Not for diagnostic use.

Please read this manual thoroughly before use.

• The antibodies used for this ELISA strongly react with various mutant viruses such as alpha, beta, gamma and delta types as well as original Wuhan type (paper is being prepared).

【 I – 3】 Kit Principle

Anti-SARS-CoV-2 spike protein antibody is immobilized on the plate. When the sample is added, the SARS-CoV-2 spike protein in the sample is trapped. After washing, the trapped protein is reacted with an HRP-labeled anti-SARS-CoV-2 spike protein antibody, and after substrate addition, the color development by HRP is read and quantified with a plate reader.

【 I – 4】 Kit Component

Storage temperature : 2 ~ 8 °C

	Reagent	Volume	Quantity
1	Anti-SARS-CoV-2 Spike protein Antibody Immobilized Plate	96well (8well x 12 strips)	1 plate
	Standard Protein	200μL (25,000pg/mL)	1 tube*
2	Assay Buffer	25mL	1 vial
3	Washing Buffer (10X)	25mL	1 vial
4	HRP Conjugated Anti- SARS-CoV-2 Spike protein Antibody (500X)	20μL	1 vial
5	Substrate Solution	12mL	1 vial
6	Stop Solution (2N H2SO4)	6mL	1 vial
8	Plate Seals		2 sheets

* Sufficient to create 4 standard curves with n=2.

Required Materials Not Included in the Kit

- Micropipettes (10 ~ 1000 μL)
- Multichannel micropipette
- Multichannel micropipette Reservoir
- Plate shaker
- Microplate reader (enable to measure at wavelength 450nm)

- Plate washer

【Ⅱ】 Preparation of Reagents and Samples

【Ⅱ – 1】 Preparation of Washing Buffer

- Dilute Washing Buffer (10×) to 10 folds with purified water.
e.g. For 1 plate, add 225mL of purified water to 25mL of Washing Buffer (10 x) and mix well.

【Ⅱ – 2】 Preparation of Standard Protein solution

	Concentration (pg/mL)	Standard Protein	Assay Buffer	Dilution factor
A	25000			
B	2500	50μL of A	450μL	10
C	1250	250μL of B	250μL	2
D	625	250μL of C	250μL	2
E	312.5	250μL of D	250μL	2
F	156.3	250μL of E	250μL	2
G	78.1	250μL of F	250μL	2
H	39.1	250μL of G	250μL	2

- To prepare Solution B, add 450μL of Assay Buffer into 50μL of Standard Protein (Solution A), and then mix well (10 times dilution). To prepare Solution C, add 250μL of Assay Buffer into 250μL of Solution B, and then mix well (2 times dilution). Similarly, 2 times dilution series for Solution D through H should be prepared.
- Use 200μL for measurement, using 2 wells for each solution (n=2).
- Diluted Standard Protein Solution(39.1~2500 ng/mL) should be freshly prepared at each time before use.

【Ⅱ – 3】 Preparation of antibody solution

- Dilute HRP conjugated anti- SARS-CoV-2 Spike protein antibody (500x) to 500 folds using Assay Buffer.

e.g. For 1 plate, add 20 μ L of antibody (500x) into 10mL of Assay Buffer. Mix by inverting the tube.

* Diluted antibody solution should be freshly prepared at each time before use.

【II –4】 Preparation of Samples

When using serum or plasma, dilute it 2-fold or more with assay buffer to make a sample solution.

【III】 Sample measurement procedure


- 1. Regarding the handling of samples and operations, take appropriate biohazard measures as necessary, and dispose of waste appropriately in accordance with each rule.**
2. Leave Anti-SARS-CoV-2 Spike protein Antibody Immobilized Plate and each reagent to the room temperature.
3. Prepare Standard Protein solution by serial dilution. (from step 【II –2】)
4. Add 100 μ L each of serial diluted Standard Protein solution (39.1 ~ 2500pg/mL) or Sample solution into the well.
5. Seal the microplate with Plate Seals, place into plate shaker, and then shake it at 800 rpm for 30 sec.
6. Incubate at room temperature for 2hrs for static reaction.
7. Discard all the reaction solution, and then rinse each well with 300 μ L of Washing Buffer (from step 【II-1】). Repeat this step for 3 times.
8. Add 100 μ L each of diluted HRP conjugated anti-SARS-CoV-2 Spike protein antibody (from step 【II-3】) to the well.
9. Seal the plate, and then shake in the plate shaker.
10. Incubate at room temperature for 2hrs for static reaction.
11. Discard the antibody solution, and then rinse each well with 300 μ L of Washing Buffer. Repeat this step for 3 times.

12. Add 100 μ L of Substrate Solution into each well, and then incubate at room temperature protected from light for 20min for static reaction.
13. Visually confirm the coloring, and then add 50 μ L each of Stop Solution.
14. Place into the plate reader, and read the absorbance of each microwell on a spectrophotometer at the wavelength of 450nm.
15. Create a standard curve for the SARS-CoV-2 Spike protein by plotting the mean absorbance (y axis) against the Standard protein concentration (x axis) (Fig. 1).
16. Calculate the concentration by comparing the absorbance obtained from the sample solution to the standard curve. Multiply the resulting value by the appropriate sample dilution factor, to obtain the concentration of SARS-CoV-2 Spike protein in the sample.

【IV】 Example of Results

【IV-1】 Standard curve

As a reference example, Fig. 1 shows the graph of absorbance (OD450) for a standard protein. However, draw a new standard curve for each assay and calculate the concentration in the sample.

Standard Protein (pg/mL)	Absorbance (450nm)		Average
	1	2	
0	0.058	0.057	0.058
39.1	0.094	0.106	0.106
78.1	0.136	0.142	0.139
156.3	0.221	0.233	0.227
312.5	0.391	0.381	0.386
625	0.738	0.697	0.718
1250	1.382	1.420	1.401
2500	2.465	2.608	2.537

Fig. 1 Standard curve and measured values

【IV-2】 Sample measurement example

Exosomes (15.625, 31.25, 62.5, 125, 250, 500, 1000ng/mL) purified from the culture supernatant of Expi293 in which SARS-CoV-2 Spike protein was forcibly expressed were measured.

Figure 2 shows the exosome concentration (ng / mL) and absorbance (OD450) graphically.


Fig. 2 Measurement of exosomes overexpressing SARS-CoV-2 spike protein (S1, S2 domain)

【Reference】

- (1) Alexandra C. Walls, Young-Jun Park, et al., Cell 180, 1-12 (2020)